



HOW TO FILL THE PRESCRIPTION AND RECOMMENDATIONS

Terms in **Black Bold** are defined in the glossary. **Purple** indicates choices recommended by Panthera Dental appliance specialists. **Red** indicates warning messages.

Patient: [Insert patient complete name]

Dentist: [Insert your name here]

Licence#: [Insert your licence practice number here]

SECTION 1

TREATMENT RANGE NEEDED (Starting Point ●)

The (●) on the pictures represents which rods your appliance will be delivered with in relation to your provided **Bite**. For each choice, 6 pairs of rods will be provided as seen on the pictures. More rods are available upon request, for an extra fee.

RETRUDE 4mm with 0.5mm step before patient's max

By selecting this option, your appliance will be optimized to start patient treatment 4mm behind your **Bite**. Other levels of protrusion can be achieved but, in rare situation, can be limited.

PROTRUDE 5mm

Select this option if you are certain to only require protrusion from your **Bite** and if you do not require retraction.

RETRUDE 1mm and Protude 4mm

Select this option if you want your appliance design to allow protrusion starting from your **Bite** but also a possibility to retrude a little from it.

VERTICAL SPACING

Close or open to optimise the device:

By selecting this option, the design team will either close or open the **Bite** to make the appliance fit and work in your patient's mouth. The team will always keep the **Bite** as closed as possible and will open only if there is no room to insert the appliance. **(The minimal Vertical Spacing required is approximately 3-4 mm)**

Keep it, call if major changes needed:

Select this option if you want your **Bite** to stay the same. In situations where the appliance cannot be inserted, the design team will open up to 2mm without any notice. If 3mm+ is needed, Panthera will contact you prior to designing the appliance.